

EIC-Therm WS 1

infusio-therm

textil-therm

Heat for patients



***Modular instrumentation system as integrated
drawers or as mobile unit for heating and warmth
retention of infusion solutions, diapers, surgical
linen or the like***

EIC-Therm WS 1

infusio-therm

textil-therm

Modular instrumentation system for heating and warmth retention of infusion solutions, diapers, surgical linen or the like.

The WS1 integrated heating drawers are compact integrated drawers in which infusion solutions, diapers, surgical linen or the like can be heated or kept in warm storage by means of infrared radiation heat, **not circulating heat**.

These integrated drawers can be installed in all existing 60 cm wide cupboard units; the front panels and handles can be designed in accordance with the model of cupboard units.

The WS1 heating drawers can be safely mounted on top of or under working platforms or trolleys with the H-WS1 universal support mechanism.

Among other things, up to twenty four 1-litre bottles or thirty six ½-litre bottles or bags of infusion solution (infusio-therm) or up to four nappy or hand towel packs (textil-therm) can be heated and stored at a temperature of up to 37°C / 55°C within 6-8 hours. The temperature can vary by +/- 2°C, depending on the type of skein and size. As a result of the infrared heating, the danger of potential germ formation through condensate water and re-circulated air heating is ruled out.

The drawers are designed as complete pull-out units.

Technical data – WS 1:

	<i>infusio-therm</i>	<i>textil-therm</i>
Voltage:	230 V	
Mains frequency:	50/60 Hz	
Nominal power:	200 W	
Conservation power:	in continuous operation: ca. 25 W at 22°C room temperature	
Fuse:	1.6 AT	
Set-point:	37°C, +/- 2°C	55°C, +/- 2°C
Temperature protection:	45°C electronic	65°C electronic
MPG-classification:	Class 1	
Protection class:	II, CE	

Drawer installation dimensions – WS 1:


Dimensions (W x H x D):	550 x 335 x 517 mm (installation dimension)
Dimensions (W x H x D):	595 x 335 x 572 mm (incl. 20 mm front panels with handle)
Filling volume (W x H x D):	480 x 265 x 408 mm
Weight:	22 kg

Variation with roller-trolley – optionally with 1 or 2 WS 1 units:

Dimension (W x H x D):	675 x 900 x 710 mm
Filling volume (W x H x D):	each WS 1, 480 x 265 x 408 mm
Rollers:	125 mm, 2 with locking lever; electrically conductive
Weight:	46 / 68 kg

Special dimensions and versions available on request

Subject to change without prior notice



Rudolf-Diesel-Straße 3
D-37308 Heilbad Heiligenstadt

Tel.: +49 3606 607925
Fax: +49 3606 5071353
E-mail: eic-ltd.mail@t-online.de